

# Steriking® COVER BAGS

**Protect Your Sterile Barrier System from Adverse Events**



## High protection

- Protects barrier and sterile item from dust, holes, abrasion, stains, insects
- Reduces level of microbes on SBS outer and inner layers
- Water repellent

## Excellent characters

- Peelable • Durable • Transparent
- Easy to seal with very wide sealing window

# Steriking® Cover bags




Steriking® Dust Cover Bags are constructed of an impermeable multilayer plastic laminate, which protects items from dust and environmental influences and thus extends the sterility maintenance time.

The STERIKING® Dust cover bags conform to the international product standards and norms: ISO 11607-1:2019 and ISO 11607-2:2019. They can be used together with wrapped trays, pouched items or containers to give extra protection for the barrier and the sterilized device. Cover bags protect sterile barrier system from breaking, high humidity and water. They provide protection against microbes, dust and dirt so that the sterility of the packed item will not be compromised.

## Cover Bags:

- CB heat sealable cover bags can be sealed with impulse or roto sealer, recommended sealing temperature 130-160° C (272-335° F)
- SSDC self-sealable cover bags are used if there is no heat sealer
- Prolongs shelf life of the sterilize items
- Cover bags have easy and clean peel opening which allows aseptic presentation



## Wipak Steriking® Cover Bags Range Includes:



STERIKING® DUST COVER BAGS HEAT SEAL			
Code	Size (mm)	Size (inch)	Packing Unit (pcs/case)
CB1	150 x 270	6' x 10,5'	1 000
CB2	205 x 300	8' x 12'	1 000
CB3	250 x 400	10' x 15,5'	1 000
CB4	300 x 500	12' x 19,5'	500
CB5	400 x 550	15,5' x 22'	500
CB6	420 x 700	16,5' x 27,5'	500
CB7	500 x 800	19,5' x 31,5'	300

STERIKING® SELF SEAL DUST COVER BAGS			
Code	Size (mm)	Size (inch)	Packing Unit (pcs/case)
SSDC10	150 x 270	6' x 10,5'	500
SSDC20	200 x 350	8' x 14'	500
SSDC30	250 x 400	10' x 15,5'	500
SSDC40	300 x 500	12' x 22'	200
SSDC60	420 x 700	16,5' x 27,5'	200

*We care that you pack safely!*



www.steriking.com  
e-mail: [steriking@wipak.com](mailto:steriking@wipak.com)

