

Flat Pouches with Tyvek®



LTS special pouches made with Tyvek® are

intended for low temperature sterilization (EO, FO, hydrogen peroxide and irradiation).

Hermetically sealable. Recommended sealing temperature 120-130° C (248-266° F).

LTS-pouches are packed in dispenser boxes of 200 pieces.

Article	Width	Length	Gusset	Packing Unit
LTS7520	75	200		1000
LTS1025	100	250		1000
LTS1530	150	300		1000
LTS1644	160	440		600
LTS1660	160	600		600
LTS2038	205	390		1200
LTS2538	250	390		600
LTS2550	250	500		600