

## Nonwoven Pouches



The Steriking®TEX-Line Peel Pouches are made of a tough, resistant nonwoven and Wipak's Multi-X9 film, offering benefits superior to peel pouches with paper packing and sheet folding process.

Added package performance:

- Stronger seals and higher mechanical strength than peel pouches with paper backing
- Increased porosity for effective air evacuation and steam or gas penetration
- Improved microbial barrier properties when compared to sheet folding techniques

More user convenience:

- Quick removal prior to delivery into use
- Safe presentation of the inner package
- Strong choice for heavier, bulkier products and sets

Improved packing process

- Less time consuming than a sheet folding process

Article	Width	Length	Gusset	Packing Unit
SNW2040	205	400		500
SNW2744	270	440		500
SNW3250	320	500		500
SNW4260	420	600		300
SNW5065	500	650		200
SNW5075	500	750		200
SNW5775	570	750		200