

Nonwoven Wraps



Nonwoven sheets with superior softness and toughness

- High tear resistance ensures safe handling of heavier items
- Use as such or in combination with crepe papers or textiles
- Economical inner packing for containers
- Good drapeability eliminates springback of sheets
- Alcohol and water resistance minimises risk of liquid-borne contamination

Nonwoven Green 52 g/m²

- Economical solution for general use, for instrument sets
- Made from virgin wood pulp and polyester fibres

Nonwoven Soft Blue 60 g/m²

- For general use, for laundry and instrument sets
- Made from virgin wood pulp and polyester fibres

Nonwoven Super Soft Green 66 g/m2

- For demanding applications
- With higher content of polyester fibres

Article	Width	Length	Gusset	Packing Unit
NWG52-60	600	600		500
NWG52-75	750	750		200
NWG52-90	900	900		200
NWG52-100	1000	1000		200
NWG52-120	1200	1200		100
NWB50	500	500		500
NWB60	600	600		500
NWB75	750	750		250
NWB90	900	900		250
NWB100	1000	1000		250
NWB120	1200	1200		100
NWB130/150	1300	1300		104
NWB137	1370	1370		100
NWG50	500	500		400
NWG60	600	600		400
NWG75	750	750		200
NWG90	900	900		150
NWG100	1000	1000		150
NWG120	1200	1200		100
NWG130/150	1300	1300		100
NWG137	1370	1370		100