

Pouches with Gussets



Gusseted pouches for more voluminous items. Printed with process indicators for steam, EO and formaldehyde sterilization.

Closing of gusseted pouches requires special attention to eliminate wrinkles in the final closure seal.

Opening of gusseted pouches calls for special care to minimize any break at points of folded edges.

Article	Width	Length	Gusset	Packing Unit
B35	100	300	50	500
B30	100	360	50	500
B31	150	400	50	500
B32	150	460	50	500
B36	200	400	55	250
B33	200	500	55	250
B37	250	480	65	250
B34	300	550	80	250
B38	420	570	100	250