

Rotosealers



ROTOSEALERS for quick and reliable closing of sterilisation packaging materials

- All-day readiness for use
- Energy saving when the machine is not in use
- All sealing process parameters are controlled
- Adjustable temperature rate
- Provide strong, impermeable seals
- Manufactured according to European safety regulations

RS120A - advanced model with alarm system

RS300 - reliable and easy to use sealing machine with integrated printer

*gusset in the product data table refers to height.

| Article | Width | Length | Gusset | Packing Unit |
|---------|-------|--------|--------|--------------|
| RS120A | 405 | 290 | 160 | |
| RS120AV | 405 | 290 | 160 | |
| RS300 | 527 | 314 | 162 | |

| Article | Width | Length | Gusset | Packing Unit |
|---------|-------|--------|--------|--------------|
| RS300V | 527 | 314 | 162 | |