

# Tray Liners



**Tray Liners for use in surgical instrument trays for absorbing moisture and for mechanical protection**

- Absorbs excessive condensation.
- Wet strengthened special paper.
- Virtually lint-free.

**TL** - standard Trayliners - into a tray under the instruments

**TLWE** - Trayliners with edge - outside of the instrument tray on the top of the inner wrap

**TLNW** - Nonwoven Trayliners - for combinations of nonwovens and SMS-type of wraps - extreme absorption capacity

Article	Width	Length	Gusset	Packing Unit
TL25	250	300		1000
TL3040	300	400		500
TL40	400	500		500

Article	Width	Length	Gusset	Packing Unit
TL50	300	500		500
TLNW1525	150	250		3200
TLNW3025	300	250		1200
TLNW3040	300	400		1200
TLNW3050	300	500		600
TLWE2525	250	250		500
TLWE4825	480	250		500
TLWE5030	500	300		500